



Nalgene® autoclavable polypropylene trays

Product Images

Description

Nalgene™ Autoclavable Polypropylene Trays Use the durable Thermo Scientific™ Nalgene™ autoclavable polypropylene trays to repeatedly autoclave instruments and laboratory equipment. Description: Autoclavable polypropylene tray Dimensions (L x W x H): 10.25 x 6.25 x 2.5 in. (260 x 159 x 64 mm) Height (Anglo): 2 1/2 in. Length (Anglo): 10 1/4 in. Quantity per package: 1 Width (Anglo): 6 1/4 in. Autoclavable: Yes Type: Autoclavable Pan Unit Size: Case of 6

Additional Information

No.	CD-1793
Manufacturer (Brand)	Nalgene

