



Nalgene® autoclavable polypropylene trays

Product Images

Description

Nalgene™ Autoclavable Polypropylene Trays Use the durable Thermo Scientific™ Nalgene™ autoclavable polypropylene trays to repeatedly autoclave instruments and laboratory equipment. Description: Autoclavable polypropylene tray Dimensions (L x W x H): 12.75 x 10.25 x 2.75 in. (324 x 260 x 70 mm) Height (Anglo): 2.75 in. Length (Anglo): 12.75 in. Quantity per Package: 1 Shape: Rectangular Width (Anglo): 10.25 in. Autoclavable: Yes Type: Pan Unit Size: Case of 6

Additional Information

No.	CD-1795
Manufacturer (Brand)	Nalgene

