



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® watch glass bowl, rim fused, d = 200 mm

€40.82 As low
as
€31.40
plus VAT &
Shipping

Product Images



Description

DURAN® watch glass dish, Ø 200 mm The multitude of different work steps in the laboratory environment demands a variety of universally applicable glassware. These glassware devices usually have a fairly simple shape and are characterized by their versatility. Watch glasses, for example, are concavely curved. The small trough created in this way allows liquids or solids to be poured in. Special features and advantages: Reduced risk of injury due to fused rim
Application: Weighing, covering, pre-sampling, drying of materials, Available in Boro 3.3 or soda-lime glass

Additional Information

No.	SD-0478
Manufacturer (Brand)	DURAN
VGKL number	112132161
Transport temperature	Room temperature

