



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## DURAN® watch glass bowl, rim fused, d = 250 mm

**€56.53**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

DURAN® watch glass dish, Ø 250 mm The multitude of different work steps in the laboratory environment demands a variety of universally applicable glassware. These glassware devices usually have a fairly simple shape and are characterized by their versatility. Watch glasses, for example, are concavely curved. The small trough thus created allows liquids or solids to be poured in. Special features and advantages: Reduced risk of injury due to fused rim  
Application: Weighing, covering, pre-sampling, drying of materials, Available in Boro 3.3 or soda-lime glass

## Additional Information

---

No.	SD-0479
Manufacturer (Brand)	DURAN
VGKL number	112132166
Transport temperature	Room temperature

