



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## DURAN® filter plate with glass rim 24 mm, NW 19, Por. 3

**€34.92**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

DURAN® filter plate, with glass rim, Ø 24 mm, porosity 3 One of the most important separation methods in the laboratory environment is filtration and is used for substance separation or purification of different media. This working method requires filters which are selected according to the application in terms of size and porosity. The filter attachment consisting of a receiver vessel, sealing ring, glass filter plate and funnel can be disassembled into its individual parts and each individual component can be replaced or varied to adapt the filter to the application. Thus, the three unit sizes combined with four filter plates each of different porosity create 12 different filter options. Special features and advantages: By using a glass filter, filter paper or membrane can be dispensed with if necessary Glass edge flamed Burned-in and therefore permanently durable white print Application: qualitative analysis and preparative chemistry, filtration without filter paper, filtration

## Additional Information

---

No.	SD-1812
Manufacturer (Brand)	DURAN
VGKL number	112520243
Transport temperature	Room temperature

