



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**DURAN® filtration apparatus acc. to Witt,  
complete with replaceable lid, Keck mounting  
set, DN 100, NS 45/40**

**€331.27  
plus VAT &  
Shipping**

## Product Images



## Description

---

DURAN® Witt filtration apparatus, sleeve NS 45/40, DN 100 Filtrations are a frequently used working method in the laboratory environment for separating a solid-liquid mixture. In filtration with the Witt filtration apparatus, the vacuum-tight vessel with its lid serves as a collecting vessel for the filtrate. A rubber sleeve is placed in the tube opening of the lid, into which the filter chute is then inserted. This working method is mainly used for collecting the filtrate. Special features and advantages: Tight seal due to standard ground joint Further ground-glass tube Burnt-in and therefore permanently durable white printed image Application: Filtration, Preferential collection of filtrate, Filtration at reduced pressure, Vacuum filtration, Safe working under vacuum

## Additional Information

---

No.	SD-1641
Manufacturer (Brand)	DURAN
VGKL number	122473146
Transport temperature	Room temperature

