



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Sartorius® Sartolab P20 Pressure Filter 18075 - D, 0.2 µm Polyethersulfone

€118.00
plus VAT &
Shipping

Product Images



Description

Sartolab® P20 is a pressure filter for sterile filtration of 100 ml to 5 L media and aqueous solutions. Main advantages:- Highest flow rates with a large filtration area (20 cm) - No protein loss with a low protein affinity membrane- Low dead volume due to optimized membrane support

Additional Information

No.	ST-3737
Manufacturer (Brand)	Sartorius
Transport temperature	Room temperature

