



Sartorius® MinisartGF,28mm,nsterile,1000pcs

Product Images



Description

Minisart plus and Minisart X plus filters contain a glass fiber prefilter in addition to CA membrane or NY(polyamide) membrane for clarification and higher throughput. Minisart® Plus syringe filters for sterile filtration of difficult-to-filter solutions. The syringe filters contain an additional glass fiber prefilter above the membrane filter. This combination has the effect of often doubling the flow rate. Minisart® Plus syringe filters are the first choice for fast sterile filtration with maximum user-friendliness. Sample purification Particle retention Sterile filtration of difficult-to-filter solutions Tested quality 100% during production

Additional Information

No.	ST-3706
Manufacturer (Brand)	Sartorius

