



Sartorius® MinisartRC, 0.45 µm, 25mm, sterile, 50pcs

Product Images



Description

Minisart RC filters with hydrophilic solvent resistant RC membrane. Ready-to-use filters for sample preparation for HPLC and GC analysis. Minisart® RC syringe filters with hydrophilic, solvent-resistant regenerated cellulose (RC) membrane. The filter units are ready-to-use and allow easy, fast and reliable ultracleaning of small sample volumes for HPLC or GC analysis. Minisart® RC syringe filters can be used for both aqueous solutions and solvent mixtures. They ensure rapid particle removal from small volumes with maximum user-friendliness and fit all standard syringes. Sample cleaning Particle retention User-friendly Tested quality 100% tested during production

Additional Information

No.	ST-3661
Manufacturer (Brand)	Sartorius

