



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bohlender® BOLA Buffer funnel NW 45 mm, 50 ml

**Price on
request
plus VAT &
Shipping**

Product Images



Description

BOLA Büchner funnel Made of PTFE. 2-part design can be disassembled into upper and lower part for easy cleaning and removal of filter cake. Suitable for vacuum filtration. Nominal sizes and dimensions according to DIN 12 905. Suitable for commercially available fi

Additional Information

No.	ER-1919
Manufacturer (Brand)	Bohlender
VGKL number	101765402
Transport temperature	Room temperature

