



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bohlender® BOLA hopper NW 90 mm, 290 ml

**Price on
request
plus VAT &
Shipping**

Product Images



Description

BOLA Büchner funnel Made of PTFE. 2-part design can be disassembled into upper and lower part for easy cleaning and removal of filter cake. Suitable for vacuum filtration. Nominal sizes and dimensions according to DIN 12 905. Suitable for commercially available fi

Additional Information

| | |
|-----------------------|------------------|
| No. | ER-1922 |
| Manufacturer (Brand) | Bohlender |
| VGKL number | 101765408 |
| Transport temperature | Room temperature |

