




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene® Two-Piece Buechner Funnel, Polypropylene

€159.00
plus VAT &
Shipping

Product Images



Description

Two-Piece Nalgene Büchner Funnels, Polypropylene Thermo Scientific Nalgene Büchner Funnels made of polypropylene are lightweight and durable. Top diameter (metric): 172 mm Capacity (metric): 2.91 l Autoclavable: Yes Description: Funnel Buchner PP 150mm Length (metric) Handle: 107 mm Quantity: 4 Unit Size: Case of 4

Additional Information

No.	CD-1489
Manufacturer (Brand)	Nalgene
Transport temperature	Room temperature

