



Nalgene® Two-Piece Buechner Funnel, Polypropylene

Product Images



Description

Two-Piece Nalgene Büchner Funnels, Polypropylene Thermo Scientific Nalgene Büchner Funnels made of polypropylene are lightweight and durable. Top diameter (metric): 130 mm Capacity (metric): 919 ml Autoclavable: Yes Description: Funnel Buchner PP 110mm Length (metric) Stem: 93 mm Quantity: 1/pk; 4/cs. Unit Size: Case of 4

Additional Information

No.	CD-1488
Manufacturer (Brand)	Nalgene

