




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene® separatory funnel, polypropylene, Teflon TFE stopcock and

€323.00
plus VAT &
Shipping

Product Images



Description

Nalgene separatory funnels, polypropylene, Teflon TFE stopcock and plug; polypropylene screw cap Thermo Scientific Nalgene separatory funnels and screw caps are chemically inert, translucent and autoclavable. Capacity (metric): 500 ml Description: 43 mm Quantity: 1/Ⓟk; 4/cs. Unit Size: Case of 4

Additional Information

No.	CD-1492
Manufacturer (Brand)	Nalgene
VGKL number	238743003
Transport temperature	Room temperature

