



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® Dropping funnel, cylindrical, with graduation, NS stopcock and stopcock safety device, 1000 ml

€185.01
plus VAT &
Shipping

Product Images



Description

DURAN® Dropping funnel, 1000 ml The metered addition of solutions, reaction mixtures and/or chemicals is ideally carried out with the aid of a dropping funnel. Special features and advantages: With NS stopcock and stopcock safety device The speed of the feed can be adjusted Cylindrical with graduation Tight closure due to standard ground joint Baked-on and therefore permanently durable white printed image Easy-to-read scale with graduation for measuring the contents Application: Uniform and metered addition of liquid to a reaction mixture, measurement of liquid quantities

Additional Information

No.	SD-1119
Manufacturer (Brand)	DURAN
Transport temperature	Room temperature

