



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® Gas wash bottle attachment with filter plate with standard ground joint 34/35, 100 ml

€225.37
plus VAT &
Shipping

Product Images



Description

DURAN® gas wash bottle, with filter plate (porosity 1), NS 34/35, 100 ml Gas wash bottles are required for gas filtration or cleaning, consisting of a gas inlet tube, a suitable bottle and a gas outlet opening. In this experimental setup, the gas is forced to flow through a liquid before it is released. The gas distribution tube is immersed in the liquid, the gas flows through it and impurities of the gas are bound in the liquid. The gas is cleaned or washed, so to speak. A filter plate in the gas inlet tube creates the smallest possible air bubbles, which thus form the largest possible contact surface. Special features and advantages: Attachment with filter plate Tight seal due to standard ground joint Burnt-in and therefore permanently durable white printed image Application: Cleaning ("washing") of gases with solvents, Cleaning of gas mixtures, The standard ground joint prevents the gas from escaping at the bond point, Filtration

Additional Information

No.	SD-1829
Manufacturer (Brand)	DURAN
VGKL number	112570101
Transport temperature	Room temperature

