



Thermo Scientific Untreated Nunc™ Flasks

Product Images

Description

Untreated Nunc Flasks Minimize the risk of contamination with these high-quality untreated Thermo Scientific Nunc flasks made of optically clear polystyrene for cell suspension cultures. Cap Color: White Description: Untreated Nunc T75 EasYFlask, Filter Cap Flask Neck Type: Angled Number Per Case: 100 Number Per Pack: 5 Cap Type: Filter Culture Area: 75 cmsup2 Surface Treatment: Untreated Unit Size: Case of 100

Additional Information

No.	FE-0165
Manufacturer (Brand)	Thermo Scientific
VGKL number	290015680
Material	Polystyrol (PS)

