



Thermo Scientific Untreated Nunc™ Flasks

Product Images

Description

Untreated Nunc Flasks Minimize the risk of contamination with these high-quality Thermo Scientific Nunc untreated flasks made of optically clear polystyrene for cell suspension cultures. Cap Color: White Description: Nunc Untreated TripleFlask, Filter Cap Flask Neck Type: Straight Number Per Case: 32 Number Per Package: 2 Cap Type: Filter Culture Area: 500 cmsup2 Surface Treatment: Untreated Unit Size: Case of 32

Additional Information

No.	FE-0032
Manufacturer (Brand)	Thermo Scientific
VGKL number	2900132903
Material	Polystyrol (PS)

