



Thermo Scientific Untreated Nunc™ Flasks

Product Images

Description

Untreated Nunc Flasks Minimize the risk of contamination with these high-quality untreated Thermo Scientific Nunc flasks made of optically clear polystyrene for cell suspension cultures. Cap Color: White Description: Untreated Nunc T175 EasYFlask, Filter Cap Flask Neck Type: Angled Number Per Case: 30 Number Per Pack: 5 Cap Type: Filter Culture Area: 175 cmsup2 Surface Treatment: Untreated Unit Size: Case of 30

Additional Information

No.	FE-0173
Manufacturer (Brand)	Thermo Scientific
Material	Polystyrol (PS)

