



Nalgene® Economy wide-mouth bottles

Product Images



Description

Nalgene Economy Wide Neck Bottles Thermo Scientific Nalgene Economy Wide Neck Bottles made of HDPE are durable and shatterproof. Capacity: 30 ml, 1 oz. Capacity (Anglo-American): 1 oz. Color: Natural, translucent Description: translucent economy wide-neck bottle Inner diameter (metric) Neck: 21 mm Number per carton: 72 Number per package: 72 Type: economy wide-neck bottle Cap color: Natural, translucent Neck finish: 28 mm Certification/Compliance: Manufactured in accordance with the quality requirements of ISO 13485:2016. Bottles and their closures meet the requirements of US standard 21CFR177.1520 for use in the food and beverage industry, USP Class VI. Closure Liner: Without liner Closure Included: Yes Autoclavable: Not autoclavable Closure Size: 28 mm Closure Material: Polypropylene Closure Type: Screw Cap Diameter (metric): 34 mm Graduated: No Height (Anglo-American): 2.49 in. Height (metric): 63 mm Opening: Wide Shape: Round Dispensing through Closure: No Width (anglo): 1.35 in. ACT Label: Yes Unit Size: Case of 72

Additional Information

No.	CD-0446
Manufacturer (Brand)	Nalgene
VGKL number	108719908

