



Nalgene® Economy wide-mouth bottles

Product Images



Description

Nalgene Economy Wide Neck Bottles Thermo Scientific Nalgene Economy Wide Neck Bottles made of HDPE are durable and shatterproof. Capacity: 500 ml, 16 oz. Capacity (Anglo-American): 16 oz. Color: Natural, translucent Description: translucent economy wide-mouth bottle Inner diameter (metric) Neck: 44 mm Quantity per carton: 48 Quantity per package: 48 Type: wide-mouth economy bottle Cap color: Natural, translucent Neck finish: 53 mm Certification/Compliance: Manufactured in accordance with the quality requirements of ISO 13485:2016. Bottles and their closures meet the requirements of US Standard 21CFR177.1520 for use in the food and beverage industry, USP Class VI. Closure Liner: Without liner Closure Included: Yes Autoclavable: Not autoclavable Closure Size: 53 mm Closure Material: Polypropylene Closure Type: Screw Cap Diameter (metric): 73 mm Graduated: No Height (Anglo): 5.52 in. Height (metric): 168 mm Opening: Wide Shape: Round Dispensing through Closure: No Width (anglo): 2.87 in. ACT Label: Yes Unit Size: Case of 48

Additional Information

No.	CD-0450
Manufacturer (Brand)	Nalgene
VGKL number	108719942

