



Nalgene® HDPE Biotainer Bottles

Product Images

Description

Nalgene HDPE Biotainer Bottles Eliminate cleaning, sterilization and related validation costs with these certified sterile, ready-to-use bottles. Capacity (metric): 4 l Capacity (anglo): 135.3 oz. Description: HDPE Biotainer, Lab Pack Diameter (anglo): 5.6 in. Inner Diameter (anglo) Neck: 1.1 in. Inner Diameter (metric) Neck: 28 mm Number Per Case: 24 Number Per Package: 8 Quantity: 8/pack 24/case Type: Standard Certification/Compliance: USP Class VI, non-cytotoxic, USP lt;661gt,, pyrogen-free, FDA 21CFR 177.1520 Closure Included: Yes Closure Size: 38 mm Diameter (metric): 143 mm For use with (devices): Storage and Transport Height (anglo): 11.77 in. Height (metric): 299 mm Opening: Narrow Shape: Round Unit Size: Case of 24

Additional Information

No.	CD-1367
Manufacturer (Brand)	Nalgene
VGKL number	108737512

