



Nalgene® IP2 wide-mouth HDPE bottles

Product Images



Description

Nalgene IP2 Wide Neck HDPE Bottles Store dry or liquid substances in Thermo Scientific™ Nalgene™ IP2 wide neck HDPE bottles, ideal as inner containers for combination packaging. Capacity: 500 ml (16 oz.) Capacity (metric): 500 ml Capacity (Anglo): 16 oz. Color:

Natural, translucent Description: wide neck bottle Passport IP2 Inner diameter (metric) Neck: 44 mm Includes: Polypropylene screw caps (PP) Height (metric) with closure: 168 mm Number per case: 48 Number per package: 12 Type: IP2 bottle Outer diameter (metric): 73 mm Certification/Compliance: Meets IATA Regulations for the Transport of Dangerous Goods, Section 5.0.2.9 Closure Included: Yes Closure Size: 53 - 415 mm Closure Material: Polypropylene (PP) Closure Type: Screw Cap Diameter (metric): 44 mm (inner), 73 mm (outer) Height (Anglo): 6.6 in. Opening: Wide Shape: Round Unit Size: Case of 48

Additional Information

No.	CD-0494
Manufacturer (Brand)	Nalgene
VGKL number	108709940

