



Nalgene® IP2 wide-mouth HDPE bottles

Product Images



Description

Nalgene IP2 wide-mouth HDPE bottles Store dry or liquid substances in Thermo Scientific™ Nalgene™ IP2 wide-mouth HDPE bottles, ideal as inner containers for combination packaging. Capacity: 1000 ml, 32 oz. Capacity (metric): 1000 ml Capacity (Anglo): 32 oz.

Color: Natural, translucent Description: wide neck bottles Passport IP2 Inner diameter (metric) Neck: 53 mm Contains: Polypropylene screw caps (PP) Height (metric) with closure: 199 mm Number per case: 24 Number per package: 6 Type: IP2 Bottle Outer diameter (metric): 91 mm Certification/Compliance: Meets IATA Regulations for the Transport of Dangerous Goods, Section 5.0.2.9 Closure Included: Yes Closure Size: 63 to 415 mm Closure Material: Polypropylene (PP) Closure Type: Screw Cap Diameter (metric): 53 mm (Inner), 91 mm (Outer) Height (Anglo): 7.8 in. Opening: Wide Shape: Round Unit Size: Case of 24

Additional Information

No.	CD-0495
Manufacturer (Brand)	Nalgene
VGKL number	108709954

