



Nalgene® Nalgene dilution bottles, polysulfone

Product Images



Description

Nalgene Dilution Bottles, Polysulfone, Polypropylene Screw Cap Thermo Scientific Nalgene dilution bottles are our most transparent and accurate dilution bottles. Capacity (metric): 220 mL Capacity (anglo): 7.37 oz. Material: polysulfone Quantity per case: 48 Quantity per package: 12 Type: Dilution Bottle Volume to Rim (anglo): 7.37 oz. Volume to Rim (metric): 220ml Closure Included: Yes Autoclavable: Autoclavable Closure Size: 38 - 415mm Closure Material: Polypropylene Closure Type: Screw Cap Graduated: No Height (Anglo): 6.05 in. Length (Anglo): 1.75 in. Opening: Wide Shape: Nalgene Square Dispensing Through Cap: No Width (Anglo): 1.75 in. Unit Size: Case of 48

Additional Information

No.	CD-0662
Manufacturer (Brand)	Nalgene

