



Nalgene® PETG Biotainer Bottles

Product Images



Description

Nalgene PETG Biotainer Bottles Save on cleaning, sterilization and related validation costs compared to reusable bottles. Capacity: 125 ml capacity (metric): 125 ml Capacity (Anglo): 4.23 oz. Description: 3025-42 Diameter (Anglo): 2 in. Inside Diameter (Anglo) Neck: 1.1 in. Inner Diameter (metric) Neck: 28 mm Number Per Carton: 100 Number Per Package: 5 Quantity: 5/pack, 70/carton Type: Square Certification/Compliance: USP Class VI, non-cytotoxic, USP It;661gt;, pyrogen-free, FDA 21CFR 177.1315 Closure Included: Yes Closure Size: 38 mm Diameter (metric): 52 mm For use with (devices): Long-term storage (application) Graduated: Yes Height (Anglo-American): 4.13 in. Height (metric): 105 mm Opening: Narrow Shape: Square Unit Size: Case of 100

Additional Information

No.	CD-0733
Manufacturer (Brand)	Nalgene
VGKL number	108730254

