



Nalgene® PETG Biotainer Bottles

Product Images



Description

Nalgene PETG Biotainer Bottles Save on cleaning, sterilization and related validation costs compared to reusable bottles. Capacity (metric): 2 l Capacity (anglo): 67.63 oz. Description: 3230-20 Diameter (anglo): 4.57 in. Inside diameter (anglo) neck: 1.4 in. Inner Diameter (metric) Neck: 36 mm Number Per Case: 20 Number Per Package: 20 Quantity: 20/pack, 20/case Type: square Certification/Compliance: USP Class VI, non-cytotoxic, USP It,661gt,, pyrogen-free, FDA 21CFR 177.1315 Closure Included: Yes Closure Size: 48 mm Diameter (metric): 116 mm For use with (devices): Long-term storage (application) Graduated: Yes Height (Anglo-American): 10.4 in. Height (metric): 265 mm Opening: Narrow Shape: Square Unit Size: Case of 20

Additional Information

No.	CD-1032
Manufacturer (Brand)	Nalgene
VGKL number	108732301

