



Nalgene® PETG Biotainer Bottles

Product Images



Description

Nalgene PETG Biotainer Bottles Save on cleaning, sterilization and related validation costs compared to reusable bottles. Capacity: 500 ml Capacity (metric): 500 ml Capacity (anglo): 16.9 oz. Diameter (Anglo): 3 in. Inner diameter (anglo) neck: 1.1 in. Inner Diameter (metric) Neck: 28 mm Number Per Package: 5 Quantity: 5/pack, 70/carton Type: square Certification/Compliance: USP Class VI, non-cytotoxic, USP It;661gt;, pyrogen-free, FDA 21CFR 177.1315 Closure Included: Yes Closure Size: 38 mm Diameter (metric): 77 mm For use with (devices): Long-term storage (application) Graduated: Yes Height (anglo): 6.93 in. Height (metric): 176 mm Opening: Narrow Shape: Square Unit Size: Case of 70

Additional Information

No.	CD-0731
Manufacturer (Brand)	Nalgene
VGKL number	108730054

