



Nalgene® PETG Biotainer Bottles

Product Images



Description

Nalgene PETG Biotainer Bottles Save on cleaning, sterilization and related validation costs compared to reusable bottles. Capacity: 5 liters (1.3 Gal.) Capacity (metric): 5 L Capacity (Anglo): 1.3 gallon. Description: With Polyethylene Handle Diameter (Anglo): 6.5 inches Inside Diameter (Anglo) Neck: 1.5 inches Inside Diameter (Metric) Neck: 38 mm Contains: Polyethylene Handle Quantity Per Carton: 6 Quantity Per Package: 6/pack, 6/carton Type: square with polyethylene handle Certification/Compliance: USP Class VI, non-cytotoxic, USP It;661gt,, pyrogen-free, FDA 21CFR 177.1315 Closure Included: Yes Closure Size: 48 mm Diameter (metric): 166 mm For use with (devices): Long-term storage (application) Graduated: Yes Height (Anglo-American): 11.4 in. Height (metric): 290 mm Opening: Narrow Shape: Square Unit Size: Case of 6

Additional Information

No.	CD-1069
Manufacturer (Brand)	Nalgene
VGKL number	108734158

