



Nalgene® Polycarbonate Biotainer bottles and

Product Images



Description

Nalgene Polycarbonate Biotainer Bottles and Balloon Vials Reduce the risk of carryover contamination and save the cost of cleaning, sterilization, and related validations. Use sterile, ready-to-use Thermo Scientific Nalgene Biotainer Bottles and Balloon Bottles, designed specifically for pharmaceutical and biological product manufacturers. Excellent for freezing. Capacity (metric): 2 l Capacity (anglo): 67.6 oz. Inside diameter (anglo) neck: 1.89 in. Inner Diameter (metric) Neck: 48 mm Quantity Per Case: 20 Type: bottle and balloon bottle Lining Material: silicone Certification/Compliance: USP Class VI, non-cytotoxic, USP lt;661 gt;, pyrogen-free, FDA 21CFR 177.1315 Closure Included: Yes Closure Size: 48 mm Closure Material: polypropylene Diameter (metric): 116 mm For use with (devices): Freezing Graduated: Yes Height (Anglo-American): 10.4 in. Height (metric): 265 mm Opening: Narrow Shape: Square Unit Size: Case of 20

Additional Information

No.	CD-1034
Manufacturer (Brand)	Nalgene
VGKL number	108732342

