



Nalgene® Polycarbonate Biotainer bottles and

Product Images



Description

Nalgene Polycarbonate Biotainer Bottles and Balloon Vials Reduce the risk of carryover contamination and save the cost of cleaning, sterilization, and related validations. Use sterile, ready-to-use Thermo Scientific Nalgene Biotainer Bottles and Balloon Bottles, designed specifically for pharmaceutical and biological product manufacturers. Excellent for freezing. Capacity (metric): 20 L Capacity (anglo): 676.3 oz. Inner Diameter (metric) Neck: 48 mm Number Per Carton: 3 Number Per Package: 3 Type: bottle and balloon bottle Liner Material: silicone Certification/Compliance: USP Class VI, non-cytotoxic, USP It;661gt,, pyrogen-free, FDA 21CFR 177.1520 Closure Included: Yes Closure Size: 48 mm Closure Material: polypropylene Diameter (metric): 255 mm For use with (devices): Freezing Graduated: Yes Height (Anglo-American): 19.4 in. Height (metric): 493 mm Opening: Narrow Shape: Square Unit Size: Case of 3

Additional Information

No.	CD-1128
Manufacturer (Brand)	Nalgene
VGKL number	108734231

