



Nalgene® Square Bottles (PP)

Product Images



Description

Nalgene Square Bottles (PP) Thermo Scientific Nalgene square bottles made of polypropylene copolymer are square and thus save space. They can be autoclaved several times. Capacity: 250 ml, 8 oz. Capacity (Anglo-American): 8 oz. Color: Natural, translucent Description: Square narrow neck bottle Diameter (Anglo-American): 2.4 in. Inside Diameter (Anglo) Neck: 1.1 in. Inside Diameter (Metric) Neck: 28 mm Includes: Polypropylene leak-proof screw caps Material: polypropylene copolymer (PPCO) Quantity per carton: 72 Quantity per package: 12 Type: Narrow Mouth Bottle Cap Color: Natural, translucent Graduations: Yes Certification/Compliance: USP Class VI; FDA - 21 CFR Closure Included: Yes Autoclavable: Yes Closure Size: 38 - 430 mm Closure Material: Polypropylene (PP) Closure Type: Screw Cap Graduated: Yes Height (Anglo-American): 5.7 in. Opening: Narrow neck Shape: Square Dispensing through closure: No Width (Anglo): 2.41 in. Unit Size: Case of 72

Additional Information

No.	CD-0289
Manufacturer (Brand)	Nalgene
VGKL number	228720164

