



Nalgene® Square Bottles (PP)

Product Images



Description

Nalgene Square Bottles (PP) Thermo Scientific Nalgene square bottles made of polypropylene copolymer are square and thus save space. They can be autoclaved several times. Capacity: 500 ml, 16 oz. Capacity (Anglo-American): 16 oz. Color: Natural, translucent Description: Square narrow neck bottle Diameter (Anglo-American): 2.9 inch inner diameter (Anglo) Neck: 1.1 inch inner diameter (metric) Neck: 28 mm Includes: Polypropylene leak-proof screw caps Material: polypropylene copolymer (PPCO) Quantity per case: 48 Quantity per package: 12 Type: Narrow Mouth Bottle Cap Color: Natural, translucent Graduations: Yes Certification/Compliance: USP Class VI; FDA - 21 CFR Closure Included: Yes Autoclavable: Yes Closure Size: 38 - 430 mm Closure Material: Polypropylene (PP) Closure Type: Screw Cap Graduated: Yes Height (Anglo-American): 6.9 in. Opening: Narrow neck Shape: Square Dispensing through closure: No Width (Anglo): 2.93 in. Unit Size: Case of 48

Additional Information

No.	CD-0290
Manufacturer (Brand)	Nalgene
VGKL number	228720165

