



Nalgene® Square Bottles (PP)

Product Images



Description

Nalgene Square Bottles (PP) Thermo Scientific Nalgene square bottles made of polypropylene copolymer are square and thus save space. They can be autoclaved multiple times. Capacity: 2 oz., 60 mL Capacity (metric): 60mL Capacity (Anglo): 2 oz. Description: Narrow-Mouth Square Bottle Diameter (Anglo-American): 1.6 in. Inside Diameter (Anglo) Neck: 0.7 in. Inside Diameter (metric) Neck: 18mm Includes: Leakproof polypropylene screw closures Quantity per case: 72 Quantity per package: 12 Type: Narrow Mouth Bottle Closure Color: Natural, translucent Graduations: Yes Certification/Compliance: USP Class VI; FDA - 21 CFR Closure Included: Yes Autoclavable: Autoclavable Closure Size: 24mm Closure Material: Polypropylene Closure Type: Screw Graduated: Yes Height (Anglo-American): 3.3 in. Opening: Narrow Shape: Square Dispensing through closure: No Width (Anglo): 1.61 in. Unit Size: Case of 72

Additional Information

No.	CD-0287
Manufacturer (Brand)	Nalgene
VGKL number	228720162

