



Nalgene® Rectangular Brown HDPE Bottles

Product Images



Description

Nalgene™ Rectangular Brown HDPE Bottles These space-saving bottles with durable walls and large opening for easy filling protect light-sensitive samples or dry materials from exposure. Capacity: 250 ml (8 oz.) Capacity (Anglo-American): 8 oz. Color: Amber Description: Brown rectangular wide-mouth bottle Inner diameter (Anglo-American) Neck: 1.1 in. Inside Diameter (metric) Neck: 28 mm Quantity Per Case: 72 Quantity Per Package: 12 Type: Wide Mouth Bottle Closure Color: Brown Certification/Compliance: Bottles meet U.S. Pharmacopeia 23 requirements for light transmission Closure Liner: Without liner Closure Included: Yes Closure Size: 38 mm Closure Material: Polypropylene (PP) Closure Type: Screw Cap Graduated: No Height (Anglo-American): 4.6 in. Opening: Wide Shape: Rectangle Dispensing through Closure: No Unit Size: Case of 72

Additional Information

No.	CD-0275
Manufacturer (Brand)	Nalgene
VGKL number	108720092

