



Nalgene® Rectangular Brown HDPE Bottles

Product Images



Description

Nalgene™ Rectangular Brown HDPE Bottles These space-saving bottles with durable walls and large opening for easy filling protect light-sensitive samples or dry materials from exposure. Capacity: 4 oz. (125 ml) capacity (Anglo-American): 4 oz. Color: Amber Description: Amber Wide-Mouth Rectangular Bottle Inside Diameter (Anglo) Neck: 0.8 in. Inner Diameter (metric) Neck: 21 mm Number Per Case: 72 Number Per Package: 12 Type: Wide Mouth Bottle Closure Color: Opaque Amber Certification/Compliance: Meet U.S. Pharmacopeia XXIII requirements for maximum light transmission. Closure Liner: Linerless Closure Included: Yes Closure Size: 28 mm Closure Material: Polypropylene Closure Type: Screw Cap Graduated: No Height (Anglo): 4 inches Opening: Wide Shape: Rectangle Dispensing through Closure: No Unit Size: Case of 72

Additional Information

No.	CD-0274
Manufacturer (Brand)	Nalgene
VGKL number	108720091

