



Nalgene® wide-mouth containers with septa

Product Images

Description

Wide-mouth jars with septa These Type III soda-lime glass wide-mouth septa jars are available in clear or amber glass and are recommended for soil analysis of volatile organic compounds. Capacity: 2 oz., 60 mL Capacity (metric): 60 mL Capacity (Anglo-American): 2 oz. Certification Parameters: VOA Description: I-Chem 60mL Amber Packer with 0.060 in. septum processed Quantity per box: 24 Replacement cap: NA Type: Jar Cap Color: White Outer Diameter (Anglo-American): 1.73 in. Fill Level: Processed Volume to Rim (Anglo): 2.4 oz. Volume to Rim (metric): 70 mL Certification/Compliance: Certified and Processed jars meet or exceed EPA performance based standards for volatile organic analysis; manufactured under an ISO certified quality management system. Closure Size: 33-400 mm Closure Material: Polypropylene with PTFE/Silicone Septum Closure Type: Open Top Cap Height (Anglo-American): 2.95 in. Unit Size: Case of 24

Additional Information

No.	CD-2220
Manufacturer (Brand)	Nalgene

