



Nalgene® wide-mouth containers with septa

Product Images

Description

Wide-mouth jars with septa These Type III soda-lime glass wide-mouth septa jars are available in clear or amber glass and are recommended for soil analysis of volatile organic compounds. Capacity: 250 ml (8 oz.) capacity (metric): 250 ml Capacity (Anglo-American): 8 oz. Certification Parameters: None Description: I-Chem 250 ml transparent tall SS jar with 0.060 inch septum, unprocessed Inner diameter (Anglo-American) Neck: 1.8 inch Inner diameter (metric) Neck: 45.7 mm Material: glass Type 3 Quantity per carton: 12 Replacement cap: S297-5800 Type: jar Cap color: White Outer diameter (Anglo-American): 2.44 inch Outer diameter (metric): 62 mm Level: Untreated Volume to rim (Anglo): 9.7 inch Volume to rim (metric): 286 ml Certification/Compliance: None Closure Size: 58 to 400 mm Closure Material: Polypropylene with PTFE/silicone septum Closure Type: Open Top Cap Height (Anglo): 5 inches Height (metric): 127 mm Opening: Wide Unit Size: Case of 12

Additional Information

No.	CD-2209
Manufacturer (Brand)	Nalgene

