



Nalgene® wide-mouth containers with septa

Product Images

Description

Wide-mouth jars with septa These Type III soda-lime glass wide-mouth septa jars are available in clear or amber glass and are recommended for soil analysis of volatile organic compounds. Capacity: 60 ml (2 oz.) Capacity (metric): 60 ml Capacity (Anglo-American): 2 oz. Certification Parameters: VOA Description: I-Chem60 ml transparent short SS jar mi t0.060 inch septum, certified Inner diameter (Anglo-American) Neck: 1.2 inch Inner diameter (metric) Neck: 30.5 mm Material: glass Type 3 Quantity per box: 24 Replacement cap: S297-5300 Type: jar Cap color: White Outer diameter (Anglo-American): 2.16 in. Outer diameter (metric): 55 mm Fill level: certified Volume to rim (Anglo): 2.5 oz. Volume to Rim (metric): 74 ml Certification/Compliance: certified and processed jars meet, at a minimum, EPA performance-based standards for volatile organic compound analysis and are manufactured in an ISO-certified quality management system. Closure Size: 53 to 400 mm Closure Material: polypropylene with PTFE/silicone septum Closure Type: open top cap Height (Anglo-American): 1.92 in. Height (metric): 49 mm Opening: Wide Unit Size: Case of 24

Additional Information

No.	CD-2228
Manufacturer (Brand)	Nalgene

