



Nalgene® wide-mouth containers with septa

Product Images

Description

Wide-mouth jars with septa These Type III soda-lime glass wide-mouth septa jars are available in clear or amber glass and are recommended for soil analysis of volatile organic compounds. Capacity: 4 oz., 120 mL capacity (metric): 120 mL capacity (anglo): 4 oz. Certification Parameters: VOA Description: I-Chem 120mL Amber short SS Jar with 0.060 in. septa certified inner diameter (Anglo-American): 1.81 in. Quantity per box: 24 Replacement cap: S297-5800 Type: Jar Cap color: White Outer diameter (anglo): 2.36 in. Fill Level: Certified Volume to Rim (Anglo): 4.6 oz. Volume to Rim (metric): 137 mL Certification/Compliance: Certified and Processed jars meet or exceed EPA performance based standards for volatile organic analysis; manufactured under an ISO certified quality management system. Closure Size: 58-400 mm Closure Material: Polypropylene with PTFE/Silicone Septum Closure Type: Open Top Cap Height (Anglo-American): 2.71 in. Opening: Wide Unit Size: Case of 24

Additional Information

No.	CD-2235
Manufacturer (Brand)	Nalgene

