



Nalgene® wide-mouth containers with septa

Product Images

Description

Wide-mouth jars with septa These Type III soda-lime glass wide-mouth septa jars are available in clear or amber glass and are recommended for soil analysis of volatile organic compounds. Capacity: 8 oz., 250 mL capacity (metric): 250 mL Capacity (Anglo-American): 8 oz. Certification Parameters: VOA Description: I-Chem 250mL Clear tall SS Jar with 0.060 in. septum certified Inside diameter (Anglo-American) Neck: 1.8 in. Inner Diameter (metric) Neck: 45.7 mm Quantity Per Case: 12 Replacement Cap: S297-5800 Type: Jar Cap Color: White Outer Diameter (Anglo): 2.44 in. Fill Level: Certified Volume to Rim (Anglo): 9.7 in. Volume to Rim (metric): 286 mL Certification/Compliance: Certified and Processed jars meet or exceed EPA performance based standards for volatile organic analysis; manufactured under an ISO certified quality management system. Closure Size: 58-400 mm Closure Material: Polypropylene with PTFE/Silicone Septum Closure Type: Open Top Cap Height (Anglo-American): 5 in. Opening: Wide Unit Size: Case of 12

Additional Information

No.	CD-2230
Manufacturer (Brand)	Nalgene

