



Nalgene® wide-neck bottles (PMP)

Product Images



Description

Nalgene Wide Neck Bottles (PMP) Thermo Scientific Nalgene clear polymethylpentene wide neck bottles can be used in the temperature range of +20 to +175C. Capacity: 8 oz., 250 mL capacity (metric): 250mL Capacity (Anglo): 8 oz. Description: Wide-Mouth Bottle Inner Diameter (metric) Neck: 33mm Quantity Per Case: 24 Quantity Per Package: 4 Type: Wide Mouth Bottle Closure Color: Natural Neck Finish: 43mm Closure Liner: Linerless Closure Included: Yes Autoclavable: Autoclavable Closure Size: 43mm Closure Material: Polypropylene Closure Type: Screw Graduated: No Height (Anglo-American): 5.05 in. Opening: Wide Shape: Round Dispensing through Closure: No Width (Anglo): 2.35 in. Unit Size: Case of 24

Additional Information

No.	CD-0369
Manufacturer (Brand)	Nalgene

