



Nalgene® wide-neck bottles (PMP)

Product Images



Description

Nalgene Wide Neck Bottles (PMP) Thermo Scientific Nalgene clear polymethylpentene wide neck bottles can be used in the temperature range of +20 to +175C. Capacity: 16 oz., 500 mL Capacity (metric): 500mL Capacity (anglo): 16 oz. Description: Wide-Mouth Bottle Inner Diameter (metric) Neck: 38mm Number Per Case: 16 Number Per Package: 4 Type: Wide Mouth Bottle Closure Color: Natural Neck Finish: 48mm Closure Liner: Linerless Closure Included: Yes Autoclavable: Autoclavable Closure Size: 48mm Closure Material: Polypropylene Closure Type: Screw Height (Anglo-American): 4.48 in. Opening: Wide Shape: Round Dispensing through Closure: No Unit Size: Case of 16

Additional Information

No.	CD-0370
Manufacturer (Brand)	Nalgene

