



Nalgene® wide-neck bottles (PMP)

Product Images



Description

Nalgene Wide Neck Bottles (PMP) Thermo Scientific Nalgene clear polymethylpentene wide neck bottles can be used in the temperature range of +20 to +175C. Capacity: 32 oz., 1000 mL capacity (metric): 1000mL Capacity (Anglo): 32 oz. Description: Wide-Mouth Bottle Inner Diameter (metric) Neck: 43mm Number Per Case: 16 Number Per Package: 4 Type: Wide Mouth Bottle Closure Color: Natural Neck Finish: 53mm Closure Liner: Linerless Closure Included: Yes Autoclavable: Autoclavable Closure Size: 53mm Closure Material: Polypropylene Closure Type: Screw Graduated: No Height (Anglo-American): 8.2 in. Opening: Wide Shape: Round Dispensing through Closure: No Width (Anglo): 3.56 in. Unit Size: Case of 16

Additional Information

No.	CD-0371
Manufacturer (Brand)	Nalgene

