



Nalgene® wide-neck bottles, square (PP)

Product Images



Description

Nalgene Wide Neck Bottles, Square (PP) Thermo Scientific Nalgene square wide neck bottles made of polypropylene copolymer feature a wide neck opening for easy cleaning, filling and removal. Capacity: 1000 ml Capacity (Anglo-American): 32 oz. Color: Natural, translucent Description: square wide-neck bottle Diameter (Anglo-American): 3.7 in. Inside diameter (anglo) neck: 2 in. Inner Diameter (metric) Neck: 53 mm Material: PPCO Quantity Per Carton: 24 Quantity Per Package: 6 Type: Square Wide Neck Bottle Closure Color: Natural Colors, Translucent Certification/Compliance: USP Class VI; FDA - 21 CFR Closure Liner: Without liner Closure Included: Yes Autoclavable: Autoclavable Closure Size: 63 mm Closure Material: Polypropylene Closure Type: Screw Cap Diameter (metric): 95 mm Graduated: No Height (Anglo): 7.1 in. Height (metric): 181 mm Opening: Wide Shape: Square Dispensing through Closure: No Width (anglo): 3.72 in. Unit Size: Case of 24

Additional Information

No.	CD-0376
Manufacturer (Brand)	Nalgene
VGKL number	108721105

