



Nalgene® wide-neck bottles, square (PP)

Product Images



Description

Nalgene Wide Neck Bottles, Square (PP) Thermo Scientific Nalgene square wide neck bottles made of polypropylene copolymer feature a wide neck opening for easy cleaning, filling and removal. Capacity: 500 ml, 16 oz. Capacity (Anglo-American): 16 oz. Color: Natural, translucent Description: square wide neck bottle Diameter (anglo): 2.9 in. inner diameter (Anglo) Neck: 1.7 in. inner diameter (metric) Neck: 43 mm Material: polypropylene copolymer (PPCO) Quantity Per Carton: 48 Quantity Per Package: 12 Type: square wide-neck bottle Closure Color: Natural, translucent Certification/Compliance: USP Class VI; FDA - 21 CFR Closure Liner: Without liner Closure Included: Yes Autoclavable: Yes Closure Size: 53 mm Closure Material: Polypropylene (PP) Closure Type: Screw Cap Graduated: No Height (Anglo-American): 5.7 in. Opening: Wide Shape: Square Dispensing through closure: No Width (Anglo): 2.97 in. Unit Size: Case of 48

Additional Information

No.	CD-0375
Manufacturer (Brand)	Nalgene
VGKL number	108721104

