



Nalgene® wide-neck bottles, square (PP)

Product Images



Description

Nalgene Wide Neck Bottles, Square (PP) Thermo Scientific Nalgene square wide neck polypropylene copolymer bottles feature a wide neck opening for easy cleaning, filling and removal. Capacity: 175 ml Capacity (Anglo): 6 oz Color: Natural, translucent Description: square wide-neck bottle Diameter (Anglo): 2 in. Inside Diameter (Anglo) Neck: 1.1 in. Inside Diameter (Metric) Neck: 28 mm Material: Polypropylene Copolymer (PPCO) Quantity Per Carton: 72 Quantity Per Package: 12 Type: Square Wide Neck Bottle Closure Color: Natural Color, Translucent Certification/Compliance: USP Class VI; FDA - 21 CFR Closure Liner: Without liner Closure Included: Yes Autoclavable: Autoclavable Closure Size: 38 mm Closure Material: Polypropylene Closure Type: Screw Cap Diameter (metric): 51 mm Graduated: No Height (Anglo): 4 inches Height (metric): 102 mm Opening: Wide Shape: Square Dispensing through closure: No Width (anglo): 2.05 in. Unit Size: Case of 72

Additional Information

No.	CD-0373
Manufacturer (Brand)	Nalgene
VGKL number	108721102

