



Nalgene® wide-neck bottles, square (PP)

Product Images



Description

Nalgene Wide Neck Bottles, Square (PP) Thermo Scientific Nalgene square wide neck bottles made of polypropylene copolymer feature a wide neck opening for easy cleaning, filling and removal. Capacity: 60 ml, 2 oz. Capacity (Anglo-American): 2 oz. Color: Natural, translucent Description: square wide-neck bottle Diameter (Anglo-American): 1.4 in. Inside Diameter (Anglo) Neck: 0.8 in. Inside Diameter (Metric) Neck: 21 mm Material: Polypropylene Copolymer (PPCO) Quantity Per Carton: 72 Quantity Per Package: 12 Type: Square Wide Neck Bottle Closure Color: Natural Color, Translucent Certification/Compliance: USP Class VI; FDA - 21 CFR Closure Liner: Without liner Closure Included: Yes Autoclavable: Autoclavable Closure Size: 28 mm Closure Material: Polypropylene Closure Type: Screw Cap Diameter (metric): 36 mm Graduated: No Height (anglo): 3.3 in. Height (metric): 83 mm Opening: Wide Shape: Square Dispensing through Closure: No Width (anglo): 1.45 in. Unit Size: Case of 72

Additional Information

No.	CD-0372
Manufacturer (Brand)	Nalgene
VGKL number	108721101

