



Nalgene® wide-neck bottles, square (PP)

Product Images



Description

Nalgene Wide Neck Bottles, Square (PP) Thermo Scientific Nalgene square wide neck bottles made of polypropylene copolymer feature a wide neck opening for easy cleaning, filling and removal. Capacity: 250 ml, 8 oz. Capacity (Anglo-American): 8 oz. Color: Natural, translucent Description: square wide-neck bottle Diameter (Anglo-American): 2.4 in. Inside diameter (anglo) neck: 1.3 in. Inner Diameter (metric) Neck: 33 mm Material: Polypropylene Copolymer (PPCO) Quantity Per Carton: 72 Quantity Per Package: 12 Type: Square Wide Neck Bottle Closure Color: Natural Colored, Translucent Certification/Compliance: USP Class VI; FDA - 21 CFR Closure Liner: Without liner Closure Included: Yes Autoclavable: Autoclavable Closure Size: 43 mm Closure Material: Polypropylene Closure Type: Screw Cap Diameter (metric): 61 mm Graduated: No Height (anglo): 4.6 in. Height (metric): 116 mm Opening: Wide Shape: Square Dispensing through Closure: No Width (anglo): 2.44 in. Unit Size: Case of 72

Additional Information

No.	CD-0374
Manufacturer (Brand)	Nalgene
VGKL number	108721103

