



Nalgene® square PETG media bottles

Product Images

Description

Nalgene square PETG media bottles Thermo Scientific Nalgene square PETG media bottles are certified sterile and non-pyrogenic. Capacity (metric): 125 mL Capacity (anglo): 4 oz. Diameter (anglo): 2.13 in. Inner diameter (Anglo) neck: 1.10 in. Inner diameter (metric) neck: 28 mm For use with (application): Sampling, Preparation, Containment, Prolonged Storage Number per carton: 50 Number per package: 5 Temperature range (metric): gt;40deg;C Certification/Compliance: USP Class VI Compliant, Certified Radiation Stabilized (SAL 10sup-6/sup) Autoclavable: Not autoclavable Closure Size: 38 to 430 mm Diameter (metric): 54 mm Height (anglo): 4.33 in. Height (metric): 110 mm Opening: Narrow Unit Size: Case of 50

Additional Information

No.	CD-3636
Manufacturer (Brand)	Nalgene

